

TEETH & GUM DISEASE CARE KIT

to CLEAN & PROTECT YOUR SMILE

SRP

DENTISTRX® REVOLUTION®

Revolving 360° Toothbrush & UV Sanitizer

- Unique power brush head design effectively reaches around tooth and braces brackets to improve the removal of plaque
- Professionally recommended home-care instrument for disease control
- Built-in timer to guide patients through the ADA recommended two-minute brushing cycle
- UV sanitizer base kills bacteria and viruses on brush head surfaces

Box Set Includes: Brush heads, easy-grip brush handle, power cord, and travel case

[More details at www.DentistRx.com]



PRO-DENRX® CHLORHEXIDINE



- Active ingredient in oral rinses to reduce dental plaque and oral bacteria
- Intended for use as an antimicrobial rinse between dental visits for the treatment of gingivitis (swelling of the gums)
- Ideal for patients undergoing active periodontal treatment

[More details at denmat.com]

PRO-DENRX® 1.1% NEUTRAL SODIUM FLUORIDE GEL



- Home-care fluoride gel
- Prescription strength - Intended for use as a topical fluoride to prevent caries (decay) and control plaque development
- Safe for use on all patients, including those with restorations

[More details at denmat.com]



westerndental.com

 **Western Dental®
& Orthodontics**

INSTRUCTIONS FOR TEETH & GUM DISEASE CARE KIT

DentistRx®



[Reference: DentistRx® Revolution® owner's manual]

DentistRx® Revolution® requires a tooth-by-tooth brushing technique rather than a scrubbing technique as described here.

- **Brushing Technique:** After applying toothpaste, hold the toothbrush at a 90-degree angle to teeth and gums, and then turn on the power. Allow the toothbrush to clean each tooth by applying light, gentle pressure for a few seconds. Then move to the next tooth and repeat. When using a revolving toothbrush, it is best to follow the gum line and use a “M” pattern when brushing your upper teeth, and a “W” pattern when brushing the lower teeth.
- **Built-In Timer:** The built-in timer will begin the 2-minute cleaning cycle pausing every 30 seconds to indicate when it's time to brush the next section of your mouth.
- **Using the UV Sanitizer:** To sanitize your brush head after use, rinse and shake dry. Place the brush head in the sanitizer tray facing the UV light and close the lid. The sanitizer will automatically run for approximately seven minutes once the lid is closed.

PRO-DENRX® CHLORHEXIDINE



Generally recommended to use for two weeks.

- First brush and floss your teeth.
- Measure out 1/2 fluid ounce (15 milliliters) as marked in the cap that comes with the bottle.
- Swish the solution in your mouth for at least 30 seconds and then spit it out.
- Do not swallow it. Do not mix it with water or other fluids. Do not eat or drink right after you use this medicine.

[More details at denmat.com]

15 milliliters
noted on inside
of bottle cap



PRO-DENRX® 1.1% NEUTRAL SODIUM FLUORIDE GEL

- Apply daily at bedtime, in place of your regular toothpaste or more often if your dentist recommends additional therapy based on the diagnosis.
- Cover brush head with Pro-DenRx 1.1% Neutral Sodium Fluoride Gel and brush around all tooth surfaces and gum line for at least two minutes.
- Spit out gel.
- **Adults:** wait 30 minutes before rinsing mouth
- **Children under 12:** Rinse mouth thoroughly immediately after use.

Cover brush
head with gel



[Reference: DenMat product user guide]



Instructions have been provided for patient convenience. Please refer to the manufacturer's product guides included in their package for details.

 **Western Dental®
& Orthodontics**